

BCS

Extraction Fume hood




Protection :
Operator

Complies with standard:
EN 14175

- ✓ Inclined sash window
- ✓ Electric glass
- ✓ Automatic Regulation ECRO Type
- ✓ Structure COPLAST
- ✓ No maintenance

The BCS is an Extraction Fume Hood with vertical flow made in COPLAST. The fume hoods ensure the protection of the user and its environment thanks to a air negative pressure of the work area.

Designed from non-particle-resistant, chemical-resistant materials, they can integrate into a controlled atmosphere zone.



	BC 9		BC 12		BC 16		BC 18	
	internal	Overall	internal	Overall	internal	Overall	internal	Overall
Width	921	957	1226	1262	1531	1567	1836	1872
Height	internal : 1200				OVERALL : 1370			
Depth	internal : 660				OVERALL : 765			