

# OPTIMALE

## Biosafety Cabinet Type II



**Triple Protection :**  
Operator  
Product  
Environment

Complies with standard:  
**EN 12469 & NSF 49**

- ✔ Electric sash window
- ✔ Automatic flow regulation ECM
- ✔ Maintenance from the front
- ✔ UV Lamp

**★ FLAGSHIP PRODUCT**

Microbiological safety stations type II are intended to protect the operator, the product and the environment from any airborne contamination.

The OPTIMALE is adapted to the handling of agents that require biosafety level class 1, 2 or 3.



	Optimale 9 NF		Optimale 12 NF		Optimale 16		Optimale 18	
	internal	Overall	Internal	Overall	internal	Overall	internal	Overall
Width	921	971	1226	1276	1531	1581	1836	1886
Height	650	1400	650	1400	650	1400	650	1400
Depth	616	786	616	786	616	786	616	786
Weight(Kg)	190		230		300		350	
Power (W)	1,4		1,4		1,5		1,5	